



DIAGNOSTIC TESTING

WHAT YOU NEED TO KNOW!

The diagnostic testing process includes three sessions.

About the sessions:

Visit 1: Initial Visit

During the first visit, the evaluator will gather a detailed history and relevant information to inform the assessment selection. This visit may take place over telehealth. Questionnaires will be sent to you after this visit and must be completed and returned before the in-person testing can be scheduled.

Additional Information

- Testing appointments are in high demand.
- Onsite testing is being scheduled in MA, NH, NY, RI, and VT. Age range and insurance accepted may vary by center. If the most convenient option falls outside of your state of residence, many insurers allow and reimburse for out-of-state diagnostics. This will require confirmation with your insurance carrier, which ACP can provide guidance on if helpful.
- A \$50 cancellation fee may be charged for no-shows or cancellations without 24 hours' notice for testing or feedback.

To Prepare For Day of Testing:

- Try to get a good night's sleep before your testing appointment.
- Continue to take all prescribed medications on the day of the appointment, unless otherwise specified by your provider, and inform your provider of all medications.
- Wear comfortable clothes and bring snacks and drinks.

Visit 2: In-Person Testing

Assessment occurs during the second visit. This can take 1-3 hours, depending on the age of the person. This visit must take place in person, and parents or guardians may be required to remain present in the center for individuals under 18.

Visit 3: Feedback

The final visit is the feedback session, where the provider will explain the results of the assessment and provide clinical impressions. Patients and their families are given the opportunity to ask any questions they may have. This visit may take place over telehealth. You will receive a copy of the diagnostic report after this appointment.

1-800-679-3609

**Diagnostic & Developmental
Evaluations Autism Care Partners**